

Bedside Lockers

Manufactured from high quality materials

A good bedside cabinet can have a positive impact on patient comfort whilst making nursing work easier. A stylish bedside cabinet is also an interior design element.

Lojer's range of bedside furniture represents quality Finnish workmanship and complements different types of beds. Bedside furniture is available in the same colour and finish options as Lojer's hospital and nursing beds.

We give our products the best warranty on the market, 2/10 years.

*"Comfortable,
stylish and high
quality."*



Reading and Meal Tray 8022

Bedside Locker 2040

Bedside Locker 2100

Bedside Locker 2020



LOJER® For easy care

■ Healthcare | Senior Care | Medical Wellness


Bed side lockers	2100	2020	2040	Reading and dining table 8022
Material, white epoxy coated steel	-	X	-	X
Material, wood	X	X	X	X
Width	40 cm	47 cm	48 cm	63 cm
Depth	46 cm	59 cm	57 cm	36 cm
Height	75 cm/ 83 cm with castors	90 cm	90 cm	67..105 cm
Weight	21 kg	35 kg	28 kg	10 kg
Height adjustment range for the tray	-	60-104 cm	60 - 104 cm	67..105 cm
Drawer dimensions, 34x37x12cm (2pcs)	X	X	X	-
Lockable wheels, 50 mm	•	X	X	-
Lockable wheels, 75 mm	•	•	•	X
Meal tray	36*42 cm	49*33 cm	49*33 cm	63*36 cm
Meal tray, tilting	-	•	•	X
Nurse calling device adapter attachment	•	•	•	-
Side basket and attachment brackets	-	•	•	-
Lock for the top drawer	•	•	•	-
Warranty 2 years	X	X	X	X
Country of Origin	Finland	Finland	Finland	Finland

Standard feature = X, Optional feature = •, Not available = -


Accessories



Meal tray, tiltable



Side basket and attachment brackets



Lock for the top drawer